



SCREENING METHOD FOR 121 ACIDIC, NEUTRAL AND BASIC DRUG ANALYTES IN PLASMA, SERUM, URINE, OR TISSUE BY LC-MS/MS

Part #:

CSXCE106 – CLEAN SCREEN XCEL[®] I 130 mg, 6 mL Tube

BETA-GLUC-10 – SELECTRAZYME[®] Beta-glucuronidase

SLDA50ID21-5UM – Selectra[®] DA HPLC Column 50 x 2.1 mm, 5 μ m

Comprehensive screening is often referred to as general unknown or systematic toxicological analysis. With the use of LC/MS/MS and full screen (QTrap methodology) the complete detection of both illicit and prescribed drugs can be accomplished easily with a simple solid phase extraction preparation. This screening tool can give low LOD's as well as specific information as to drugs and metabolites in an unknown sample. The objective of this screening application is to develop a rapid, highly sensitive qualitative method for the analysis of acidic, neutral, and basic compounds in biological fluids using Clean Screen XCEL[®] I solid phase extraction columns. The sample preparation is minimized while efficiently extracting a large group of representative compounds. Sample analysis was executed using a unique polyaromatic phase HPLC column.

1. PREPARE SAMPLE:

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards

Add 1 -2 mL of blood, plasma/ serum, urine, or 1 g (1:4) tissue homogenate

Mix/vortex and let stand for 5 minutes

Add 2 mL of 100 mM phosphate buffer (pH 6.0). Mix/vortex

Sample pH should be 6.0 \pm 0.5.

Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate. Centrifuge for 10 minutes at 2000 rpm and discard pellet

NOTE: See Hydrolysis step if required

Hydrolysis (for urine samples only): To 1-2 mL of urine sample, add 1 mL of acetate buffer

(pH 5.0) containing 5,000 units/mL Selectrazyme[®] β -glucuronidase.

Optionally, add 1 mL of acetate buffer and 25-50 μ L of concentrated β -glucuronidase.

Vortex and heat for 1-2 hours at 65°C.

2. APPLY SAMPLE:

Load sample directly to column without any preconditioning

Pull sample through at a rate of 1-2 mL/ minute

Dry column thoroughly under full vacuum or positive pressure for 1 minute

3. WASH 1 – ACIDIC & NEUTRAL COMPOUNDS (FRACTION 1):

Add 1 x 1 mL of DI H₂O

Apply pressure to column for ~1minute (either vacuum (10mm Hg) or positive pressure(~80-100psi). This ensures that the entire sample and any residual is pulled through to waste

Add 1 x 1 mL of 0.1M Acetic Acid

Apply pressure to column for ~1minute (either vacuum (10mm Hg) or positive pressure(~80-100psi).

Add 1 x 2 mL Hexane to remove residual aqueous phase

Dry column (5 minutes at full vacuum or pressure)

4. ELUTION 1 – ACIDIC & NEUTRAL COMPOUNDS (FRACTION 1):

Add 1 x 1 mL Ethyl Acetate: Hexane (50:50) Collect eluate at 1 to 2 mL/minute

5. DRY ELUTE:

Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C

Reconstitute with 100 µL of Ethyl Acetate or Mobile Phase

6. WASH 2 - BASIC COMPOUNDS (FRACTION 2):

Add 1 x 1 mL of 2% Acetic Acid/98% Methanol

Dry column 5 minutes at full vacuum (10mm Hg) or positive pressure (~80-100 psi)

7. ELUTION 2 - BASIC COMPOUNDS (FRACTION 2):

1 x 1 mL of CH₂Cl₂/ IPA/ Ammonium Hydroxide (78/20/2).

8. DRY ELUTE:

Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C. Take care not to overheat or over evaporate. Certain compounds are heat labile, such as the amphetamines. Reconstitute with 100 µL Mobile Phase.

Notes:

(1) Fraction 1 (Acid Neutrals) and Fraction 2 (Bases) can be combined together if need be. This is not generally recommended as the Acid/ Neutral fraction tends to be dirtier than the Basic one, so for more effective results, keep fractions separate.

(2) A keeper solvent such as DMF can be used to prevent the volatilization of amphetamines. Use 30-50 µL of high purity DMF in the sample (Fraction 2) before evaporation.

(3) A 1% HCl in CH₃OH solution has been used to prevent volatilization by the formation of the hydrochloric salt of the drugs. Add 1 drop of the solution prior to evaporating then continue to dryness.

(4) The hexane wash step can be removed if user is looking to analyze for Parent THC.

(5) To extract the benzodiazepine group at higher recovery, following the elution of the Acidic/Neutral drugs, a second elution can be done prior to the second wash phase. The second elution solvent consists of 98% Ethyl Acetate/ 2% Ammonium Hydroxide.

INSTRUMENT CONDITIONS (LC-MS/MS):

PARAMETERS

Instrument: Shimadzu HPLC 20-AD

Detector: AB Sciex API 3200 Qtrap MS/MS

LC Column: UCT Selectra[®] DA HPLC Column 50 x 2.1mm, 5 µm

Polarity: ESI +

Mobile Phase A: 0.1% Formic Acid in D.I. H₂O

Mobile Phase B: 0.1% Formic Acid in Methanol


Flow Rate: 0.6 mL/minute

Injection Volume: 20 µL

Gradient:

Time	%A	%B
0.00	90	10
0.50	90	100
4.00	60	40
7.50	15	85
8.50	10	90
8.51	90	10
10.00	STOP	

CHROMATOGRAM



ANALYTE TABLE

ANALYTE	Relative Retention Time (min)	Q1	Q3	LOD (ng/mL) (blood samples)
Ecgoninemethylester	0.5	200.1	182.1	5
Phenylpropanolamine	0.9	152.2	134.2	20
Morphine	1.4	286	152	2
Oxymorphone	1.5	302	227	2
Pregabalin	1.5	160.2	97	50
Pseudoephedrine	1.9	166.1	148.1	20
Hydromorphone	1.9	286	185	2
Ephedrine	1.9	166.2	148.3	20
Amphetamine	2	136.1	91.1	10
Acetaminophen	2	152	110	50
Gabapentin	2.2	172.1	67.1	50
3,4-Methylenedioxyamphetamine	2.5	180.1	105	5
Atropine	2.5	290.2	124.1	50
Buspirone	2.5	386.2	122.1	20
Clonidine	2.5	230	213	10
Methamphetamine	2.5	150.1	91.1	10
Nicotine	2.5	163.1	132.1	20
Phenylephrine	2.5	168.1	91.1	20
Theobromine	2.5	181.1	138	20
Theophylline	2.5	181.1	124	20
Mephedrone	2.5	178.2	160.1	10
Phentermine	2.5	150.2	91.2	10
6-O-Monoacetylmorphine	2.6	328.1	165.1	2
Naloxone	2.8	328.2	310.2	2
Methylone	2.8	208	160.1	10
Phenmetrazine	2.8	178.2	115.1	20
Phendimetrazine	2.8	192.2	147.1	20
Caffeine	3	195.1	122.9	50
Dihydrocodeine	3	302.2	199.1	2
Codeine	3	300	152	2
Desmethyltramadol	3	250.2	58.2	5
3,4-Methylenedioxymethamphetamine	3.1	194.1	105.1	10
7-Aminonitrazepam	3.1	252.1	121.1	10
Oxycodone D6	3.1	322.3	304.1	NA
Oxycodone	3.2	316.1	298.1	2
Hydrocodone	3.4	300	199	2
Diethylpropion	3.4	206.2	100.2	20
3,4-Methylenedioxyethylamphetamine	3.6	208.1	77.1	10
Naltrexol	3.6	344.3	308.4	5
Pheniramine	3.8	241.2	167.2	10
Olanzapine	4	313.1	256.1	20

ANALYTE	Relative Retention Time (min)	Q1	Q3	LOD (ng/mL)
Norketamine	4	224.1	207.1	10
Methylphenidate	4.1	234.1	84.1	20
Norfentanyl	4.1	233.2	84.1	5
Doxylamine	4.1	271.3	167.2	20
Nalbuphine	4.1	358.4	185.2	5
Tramadol	4.3	264.2	58	5
Tapentadol	4.3	222.3	107.2	20
Benzoylcegonine	4.4	290.1	168.1	5
7-Aminoclonazepam	4.5	286.1	121.1	5
Ketamine	4.5	238.1	125	10
Meperidine	4.5	248.2	220	20
Meprobamate	4.6	219.1	158.2	25
Normeperidine	4.7	234.1	91.2	20
Cocaine	4.9	304.1	182.1	5
MDPV	5	276.2	126.2	10
Midazolam	5	326.1	291.3	10
Bupropion	5	240.2	184	20
alpha-pyrrolidinopentophenone	5	272.3	110.1	10
5-methoxy DALT	5	272.3	110	10
7-Aminoflunitrazepam	5.2	284.1	135.1	10
Chlorpheniramine	5.2	275.1	230.1	20
Venlafaxine	5.2	278.2	260.2	25
Mirtazapine	5.3	266.2	195.1	10
Pentazocine	5.3	286.3	175.1	5
Norbuprenorphine	5.4	414.2	187.1	5
Butorphanol	5.4	328.4	131.2	5
Brompheniramine	5.5	319.1	274.1	20
Clozapine	5.5	327.1	270.1	20
Zolpidem	5.6	308.2	235.2	20
Diphenhydramine	5.8	256.2	165.1	10
Buprenorphine	5.8	468.2	396.2	5
Citalopram	5.9	325.2	109	10
D3-Doxepin	5.9	283	107.1	NA
Trazodone	5.9	372.2	176.1	5
Doxepin	6	280.2	107.1	10
Fentanyl	6	337.2	188.2	1
Fluoxetine	6	310.1	117.1	20
Haloperidol	6	376.1	123	10
Clomipramine	6	315.2	86.1	10
Phencyclidine-D5	6	249.2	164.2	NA
Dextromethorphan	6.1	272.2	171.2	5
Mianserin	6.1	265.2	208.2	20
Phencyclidine	6.1	244.2	86.1	5
Carisoprodol	6.1	261.2	176.1	100

ANALYTE	Relative Retention Time (min)	Q1	Q3	LOD (ng/mL)
Quetiapine	6.2	384.2	253.1	20
Zopiclone	6.2	389.1	245	25
Dextropropoxyphene	6.3	340.2	266.2	15
Propoxyphene	6.3	340	58	10
alpha-hydroxymidazolam	6.3	342.1	168.1	5
Desipramine	6.4	267.2	72.1	10
Imipramine	6.4	281.2	86.1	20
EDDP	6.4	278.2	234.1	20
Cyclobenzaprine	6.4	276.2	215	10
Bromazepam	6.5	316	182.1	20
Nortriptyline	6.5	264.2	233.1	20
Paroxetine	6.5	330.1	192.1	50
Carbamazepine	6.5	237.1	194.2	50
Amitriptyline	6.6	278.2	233.2	10
Lorazepam	6.8	321	229.1	10
Methadone	6.8	310.2	265.2	5
Clonazepam	6.9	316.1	270.1	10
Desalkylflurazepam	6.9	289	140.1	10
Oxazepam	6.9	287.1	241.1	5
alpha-Hydroxytriazolam	6.9	359	331.1	10
2-Hydroxyethylflurazepam	7	333.1	211.2	10
Chlordiazepoxide	7	300.1	227.1	10
Triazolam	7	343	239	10
alpha-Hydroxyalprazolam	7	325.1	297.2	5
Norfluoxetine	7	296.2	134.2	50
Nordiazepam	7.2	271.1	140.1	5
Sertraline	7.2	306.1	159	20
Estazolam	7.3	295.1	205.2	5
Flunitrazepam	7.3	314.1	268.1	5
Alprazolam-D5	7.3	314.2	286.3	NA
Alprazolam	7.4	309.1	281.1	5
Temazepam	7.4	301.1	255.1	10
D5-Diazepam	7.5	290	198.2	NA
Diazepam	7.7	285.1	193.2	5
Methaqualone-d7	8	259.2	98.2	NA
Flurazepam	8.3	388.1	315.1	5

5102-04-04