



## Analysis of Cyromazine in Poultry Feed Using a QuEChERS Approach

UCT Product Number:

**ECMSSA50CT-MP** (6 g anhydrous MgSO<sub>4</sub> and 1.5 g Na Acetate)

**EEC18156** (500 mg endcapped C18, 6 mL cartridge)

October 2010

### Introduction

This summary outlines a QuEChERS procedure for the analysis of the insecticide cyromazine (Trigard or Larvadex) in poultry feed by LC-MS/MS. Processing time is significantly faster than EPA method AG-555 and uses less solvent. Modifications include adding glacial acetic acid to the acetonitrile to increase extraction efficiency.

### Procedure

#### 1. Sample Preparation

- a) Homogenize 2 grams of poultry feed and add to a 50 mL centrifuge tube
- b) Add 10 mL of acetonitrile/acetic acid (75:25)
- c) Sonicate at 50/60 Hz for 15 minutes
- d) Add the contents of **ECMSSA50CT-MP** pouch and shake for 1 minute
- e) Centrifuge at 3400 rpm for 10 minutes
- f) Transfer 1 mL of supernatant to a calibrated test tube and add 9 mL of water: acetonitrile (95:5) with 0.1% acetic acid

#### 2. Sample Clean-up

- a) Add the 10 mL from 1) f) above to a **EEC18156** cartridge and elute dropwise
- b) Filter eluant using a 0.45 µm Teflon filter (Millipore, Billerica, MA) or equivalent
- c) Transfer 2 mL of eluant to an HPLC vial for analysis by LC-MS/MS

#### 3. Analysis LC-MS/MS

Waters Alliance 2695 HPLC (Waters) coupled with a micromass Quattro Micro triple-quadrupole mass spectrometer (Micromass, Manchester, U.K.) or equivalent

**HPLC conditions:****Guard column** (Alltima, C18, 5 µm, 2.1 x 7.5 mm, Deerfield, IL) or equivalent**Analytical column** (Alltima, C18, 5 µm, 2.1 x 250 mm, Waters) or equivalent**Mobile phase:** (A) acetonitrile with 0.1% formic acid and (B) water with 0.1% formic acid**Gradient:**

- 0-2 min, 5%A
- 2-5 min from 5 to 10% A
- 5-5.5 min from 10 to 90% A
- 5.5-8 min 90 to 5% A
- 8-10 min, from 90 to 5% A
- 10-12 min, 5% A

**Flow rate** 0.2 mL/minute**Injection volume:** 25 µL**• Mass Spectrometer**

- Positive ion mode electrospray ionization
- Monitor the ion transition of the parent ion (m/z 167) to the product ion (m/z 85) in multiple reaction monitoring (MRM)

**Mass Spectrometry Conditions for Cyromazine Quantitation**

<b>Capillary Voltage</b>	3.1 kV
<b>Cone Voltage</b>	65 V
<b>Collision Energy</b>	21-24 V
<b>Source Temperature</b>	120° C
<b>Desolvation Temperature</b>	350° C
<b>Cone Gas Flow</b>	135 L/h
<b>Desolvation Gas Flow Rate</b>	750 L/h
<b>Collision Gas</b>	Argon
<b>Parent Ion</b>	(m/z) 167
<b>Product Ion</b>	(m/z) 85

\*Summarized with permission from Xia, Kang, Atkins, Jack et al, "Analysis of Cyromazine in Poultry Feed Using the QuEChERS Method Coupled with LC-MS/MS" J. Agric. Food Chem, DOI:10.1021/jf9034282

Listing of instrument manufacturers does not constitute endorsement by UCT

DCN-012101-196