



Safety Data Sheet

Enviro-Clean 08 Endcapped 500mg/3mL

Section 1. Identification

Product Identifier Enviro-Clean 08 Endcapped 500mg/3mL
Synonyms N/A
Manufacturer Stock Numbers EEC08153

Recommended use N/A
Uses advised against N/A

Manufacturer Contact Address
UCT LLC
2731 Bartram Road
Bristol, PA, 19007
US

Phone
(215) 781-9255 0

Emergency Phone Fax
(800) 424-9300 Chemtrec N/A
CHEMTREC

Email
manastasio@unitedchem.com

Website
www.unitedchem.com

Section 2. Hazards Identification

Classification N/A
Signal Word Warning
Pictogram



Hazard Statements
Causes eye irritation
May be harmful if swallowed
May cause an allergic skin reaction
May cause respiratory irritation

Precautionary Statements

Response	If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If inhaled: Remove person to fresh air and keep comfortable for breathing. If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If swallowed: Call a poison center/doctor if you feel unwell.
Prevention	N/A
Storage	N/A
Disposal	N/A
General	If medical advice is needed, have product container or label at hand Read label before use
Ingredients of unknown toxicity	0%
Hazards not Otherwise Classified	
Medical Conditions Aggravated	Perclude from exposure those individuals having a history of respiratory illness or pre-exisiting eye or skin conditions.
Additional Information	May cause irritation to the eyes, skin and respiratory tract.

Section 3. Ingredients

CAS	Ingredient Name	Weight %
	Organically modified bonded Silica	Unknown

Occupational exposure limits, if available, are listed in Section 8.

Section 4. First-Aid Measures

General	Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area
Eye	Immediately flush eyes with plenty of water. Get medical attention, if irritation persists.
Skin	Scrub with soap and water. Remove contaminated clothing and shoes.
Ingestion	Get medical attention
Inhalation	Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention.
Comments	If symptoms persist, get medical attention. Never give anything by mouth to an unconscious person.

Section 5. Fire Fighting Measures

Suitable Extinguishing Media	Not considered to be a flammable hazard. Not considered to be an explosion hazard. water spray, CO2, foam, dry chemical As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.
Unsuitable Extinguishing Media	N/A

Section 6. Accidental Release Measures

General Procedures	Ventilate area of leak or spill. Wear appropriate personal protective equipment as specified in Section 8. Spills: Sweep up and containerize for reclamation or disposal. Vacuuming or wet sweeping may be used to avoid dust dispersal.
Special Protective Equipment	Wear appropriate personal protective equipment as specified in section 8.

Section 7. Handling and Storage

Handling	Use with adequate ventilation.
Storage.	Store in a cool and dry place. Store in a tightly sealed container. Protect from heat, direct sunlight and source of ignition. Store away from acidic, alkaline and oxidizing materials. Provide adequate ventilation.
Additional Information	If pouring into a flammable liquid, ground both containers

Section 8. Exposure Controls/Personal Protection

Occupational Exposure Limits	Ingredient Name	ACGIH TLV	OSHA PEL	STEL
	Organically modified bonded Silica	N/A	N/A	N/A
Personal Protective Equipment	Goggles, Gloves, Respirator			
Engineering Controls	Facilities storing or utilizing this material should be equipped with an eyewash station and a safety shower. Use process enclosures, local exhaust ventilation, or other engineering controls to control airborne levels below recommended exposure limits.			
Eye and Face	Use chemical goggles			
Skin	Use impervious gloves. Use impervious clothing as necessary to protect against skin contact			
Respiratory	Maintain adequate ventilation. Mechanical ventilation recommended. In case of exposure, use appropriate NIOSH approved respiratory protection.			
Work Hygienic Practices	General industrial hygiene practice			
Comments	To identify additional Personal Protective Equipment (PPE) requirements, it is recommended that a hazard assessment, in accordance with the OSHA PPE Standard (29 CFR 1910.132), be conducted before using this product.			

Section 9. Physical and Chemical Properties

Physical State	Solid
Color	White powder
Odor	N/A
Odor Threshold	N/A
Solubility	N/A
Partition coefficient Water/n-octanol	N/A
VOC%	N/A
Viscosity	N/A
Specific Gravity	1
Density lbs/Gal	N/A
Pounds per Cubic Foot	N/A
Flash Point	N/A
FP Method	N/A
Ph	N/A
Melting Point	N/A
Boiling Point	N/A
Boiling Range	N/A
LEL	N/A
UEL	N/A
Evaporation Rate	N/A
Flammability	N/A
Decomposition Temperature	N/A
Auto-ignition Temperature	N/A
Vapor Pressure	N/A
Vapor Density	N/A

Physical State Comments All data is based on bonded silica gel. Bonded silica gels have not been completely tested for health effects. This product contains synthetic amorphous silica, not to be confused with crystalline silica such as quartz, cristobalite or tridymite or with diatomaceous earth or other naturally occurring forms of amorphous silica that frequently contain crystalline forms.

Section 10. Stability and Reactivity

Stability (Under Normal Conditions) Stable under ordinary conditions of use and storage
Polymerization Hazardous polymerization will not occur
Conditions to Avoid Heat, moisture and incompatibles
Incompatible Materials Reacts with Hydrogen flouride, flourine, Oxygen diflouride, Chlorine triflouride, strong acids, strong bases and oxidizers.

Section 11. Toxicological Information

Comments No toxicity information available

Section 12. Ecological Information

Environmental Data	When released into the soil, this material is not expected to biodegrade
Aquatic Toxicity	This material is not expected to be toxic to aquatic life

Section 13. Disposal

Method	Incinerate. Follow all federal, state and local regulations
Comments	<p>It is recommended that waste generators determine whether a discarded chemical is classified as a hazardous waste.</p> <p>Additionally, waste generators must consult with their waste broker and state and local hazardous waste regulations to insure complete and accurate classification.</p> <p>Although not listed RCRA hazardous waste, this material may exhibit one or more characteristics of a hazardous waste and require appropriate analysis to determine specific disposal requirements.</p> <p>US EPA guidelines for classification determination are listed on 40CFR Parts 261.3</p>
Empty Container	Containers of this material may be hazardous when empty since they retain product residues (liquids). Observe all warnings and precautions listed for this product.

Section 14. Transport Information

UN Number	N/A
UN Proper Shipping Name	Not a Dangerous Goods
DOT Classification	N/A
Packing Group	N/A
Transport	Not Dangerous Goods

Section 15. Regulatory Information

TSCA	<p>United States</p> <p>TSCA (Toxic Substance Control Act)</p> <p>Status: This product is listed on the TSCA inventory</p>
Additional Information	<p>Canada</p> <p>Domestic Substance List (Inventory): CAS #71889-02-6 is listed on Canada's DSL list.</p>

Section 16. Other Information

Revision Date

2/22/2016

Manufacturer Supplemental Notes

All technical data is believed to be correct but does not purport to be all inclusive and shall be used only as a guide.

Manufacturer Disclaimer

The data in this Safety Data Sheet relates only to the specific material designated herein. It does not relate to use in combination with any other material or in any process.

Statement

THE INFORMATION AND RECOMMENDATIONS CONTAINED IN THIS SAFETY DATA SHEET HAVE BEEN COMPILED FROM SOURCES BELIEVED TO BE RELIABLE AND TO REPRESENT THE BEST OPINION ON THE SUBJECT AS OF THE DATE ON THIS SHEET. HOWEVER, NO WARRANTY, GUARENTEE OR REPRESENTATION, EXPRESSED OR IMPLIED, IS MADE BY UCT, INC. AS TO THE CORRECTNESS OR SUFFICIENCY OF THIS INFORMATION OR TO THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.

Additional Information

HMIS Rating
Health - 2
Flammability - 1
Physical Hazard - 1

NFPA Codes
Health Hazard - 1
FLammability - 0
Reactivity - 0