




# Spin Filtr<sup>TM</sup>

A Whole **New Spin** on dSPE





## UCT's new SpinFiltr™ is taking the hassle out of dSPE

- Simply discard dSPE chamber containing unwanted matrix and sorbent following centrifugation.
- Recover additional sample volume without worrying about the disruption of centrifugation layers.
- Obtain enhanced extract purification with our built-in 0.2 µm filtration device.

## Discover a whole new approach to dSPE!

Product Number	Contents
ECQUSF14CT	150 mg MgSO <sub>4</sub> 50 mg PSA
ECQUSF24CT	150 mg MgSO <sub>4</sub> 50 mg PSA 50 mg C18
ECQUSF34CT	150 mg MgSO <sub>4</sub> 25 mg PSA 2.5 mg GCB
ECQUSF44CT	150 mg MgSO <sub>4</sub> 50 mg PSA 50 mg C18 7.5 mg GCB
ECQUSF54CT	150 mg MgSO <sub>4</sub> 50 mg PSA 50 mg C18 50 mg Chlorofiltr®
ECQUSF64CT	150 mg MgSO <sub>4</sub> 50 mg C18 50 mg Chlorofiltr®

